


# Polymerase Chain Reaction (PCR)



# PCR: Amplification of a DNA fragment




- Portion of chloroplast DNA
- Portion of mitochondrial DNA
- Portion of bacterial DNA, etc.

# PCR Equipment and Materials

Vortex



Racks



Pipettes



p200



p20



Pipette tips

PCR tube strips



Marker



PCR machine



PCR reagents



# PCR reagents and PCR protocol

PCR reagents:  
PCR PreMix  
Primers 1 and 2  
water

Follow the protocol

## PCR protocol

### (plant DNA detection from leafhoppers)

1. Turn the thermocycler on (the black one on the left)
2. Take out from the freezer: PCR PreMix, primers (working solution; 4  $\mu$ M)
3. Prepare on the bench: ddH<sub>2</sub>O (in the right drawer), one 1.5 ml microcentrifuge tube.
4. Prepare PCR cocktail (#reactions = #samples+1; to cover pipetting error):

PCR reagent	1 reaction	5 reactions	10 reactions
PCR PreMix	10 $\mu$ l	50 $\mu$ l	100 $\mu$ l
Primer 1 (2 $\mu$ M)	2 $\mu$ l	10 $\mu$ l	20 $\mu$ l
Primer 2 (2 $\mu$ M)	2 $\mu$ l	10 $\mu$ l	20 $\mu$ l
ddH <sub>2</sub> O	5.2 $\mu$ l	26 $\mu$ l	52 $\mu$ l

5. Prepare 0.2 ml PCR strip tubes (in the drawer, next to the window)
6. Take DNA samples out from the fridge.
7. Place **19.2  $\mu$ l** of PCR cocktail in each tube.
8. Add **0.8  $\mu$ l** of a DNA sample in each tube. Mixed by pipetting.
9. Place strip tubes in the thermocycler (in any block: right or left)
10. Click on "RUN" -> "PROCEED". "Under "MAIN" find protocol ALINA1LH. Run it.

#### PCR settings:


- 94°C for 4 min
- $\times 35$  cycles:
  - 94°C for 30 sec
  - 57°C for 30 sec
  - 72°C for 30 sec
- 72°C for 2 min
- Hold at 4°C for 9 hours

← Typically, PCR runs about 2-2.5 hours


11. To stop the protocol, press on "CANCEL". Turn the thermocycler off. Take the tubes out and place them in the freezer.

# Preparing PCR “cocktail”

PCR PreMix



DNase/RNase-  
Free Distilled  
Water



Primer 1 (forward)

Primer 2 (reverse)



10  $\mu$ l

5.2  $\mu$ l

2  $\mu$ l

2  $\mu$ l



## PCR protocol

(plant DNA detection from leafhoppers)

1. Turn the thermocycler on (the black one on the left)
2. Take out from the freezer: PCR PreMix, primers (working solution: 4  $\mu$ M)
3. Prepare on the bench: ddH<sub>2</sub>O (in the right drawer), one 1.5 ml microcentrifuge tube.
4. Prepare PCR cocktail (#reactions = #samples+1; to cover pipetting error):

PCR reagent	1 reaction	5 reactions	10 reactions
PCR PreMix	10 $\mu$ l	50 $\mu$ l	100 $\mu$ l
Primer 1 (2 $\mu$ M)	2 $\mu$ l	10 $\mu$ l	20 $\mu$ l
Primer 2 (2 $\mu$ M)	2 $\mu$ l	10 $\mu$ l	20 $\mu$ l
ddH <sub>2</sub> O	5.2 $\mu$ l	26 $\mu$ l	52 $\mu$ l

# Demo 1: Adding PCR PreMix and Water

# Demo 2: Adding Primers



Alina Avanesyan

# Before the next step...

After completing step 4, PCR reagents should be placed back to the freezer



## PCR protocol

### (plant DNA detection from leafhoppers)

1. Turn the thermocycler on (the black one on the left)
2. Take out from the freezer: PCR PreMix, primers (working solution; 4  $\mu$ M)
3. Prepare on the bench: ddH<sub>2</sub>O (in the right drawer), one 1.5 ml microcentrifuge tube.
4. Prepare PCR cocktail (#reactions = #samples+1; to cover pipetting error):

PCR reagent	1 reaction	5 reactions	10 reactions
PCR PreMix	10 $\mu$ l	50 $\mu$ l	100 $\mu$ l
Primer 1 (2 $\mu$ M)	2 $\mu$ l	10 $\mu$ l	20 $\mu$ l
Primer 2 (2 $\mu$ M)	2 $\mu$ l	10 $\mu$ l	20 $\mu$ l
ddH <sub>2</sub> O	5.2 $\mu$ l	26 $\mu$ l	52 $\mu$ l

5. Prepare 0.2 ml PCR strip tubes (in the drawer, next to the window)
6. Take DNA samples out from the fridge.
7. Place **19.2  $\mu$ l** of PCR cocktail in each tube.
8. Add **0.8  $\mu$ l** of a DNA sample in each tube. Mixed by pipetting.
9. Place strip tubes in the thermocycler (in any block: right or left)
10. Click on "RUN" -> "PROCEED". "Under "MAIN" find protocol ALINA1LH. Run it.

#### PCR settings:


- 94°C for 4 min
- $\times 35$  cycles:
  - 94°C for 30 sec
  - 57°C for 30 sec
  - 72°C for 30 sec
- 72°C for 2 min
- Hold at 4°C for 9 hours

11. To stop the protocol, press on "CANCEL". Turn the thermocycler off. Take the tubes out and place them in the freezer.



# Adding DNA Template

PCR PreMix



DNase/RNase-  
Free Distilled  
Water



Primer 1 (forward)

Primer 2 (reverse)




10  $\mu$ l


5.2  $\mu$ l

2  $\mu$ l

2  $\mu$ l



19.2  $\mu$ l (in each tube)



0.8  $\mu$ l




5. Prepare 0.2 ml PCR strip tubes (in the drawer, next to the window)
6. Take DNA samples out from the fridge.
7. Place **19.2  $\mu$ l** of PCR cocktail in each tube.
8. Add **0.8  $\mu$ l** of a DNA sample in each tube. Mixed by pipetting.

# Demo 3: Transferring PCR “cocktail” to PCR tubes



Alina Avanesyan

# Demo 4: Adding DNA Template





Alina Avanesyan

# Demo 5: Labelling PCR tubes



Alina Avanesyan

# Running PCR



~2 hours later...



9. Place strip tubes in the thermocycler (in any block: right or left)
10. Click on "RUN" -> "PROCEED". "Under "MAIN" find protocol ALINA1LH. Run it.

PCR settings:

- 94°C for 4 min
- ×35 cycles:
  - 94°C for 30 sec
  - 57°C for 30 sec
  - 72°C for 30 sec
- 72°C for 2 min
- Hold at 4°C for 9 hours

11. To stop the protocol, press on "CANCEL". Turn the thermocycler off. Take the tubes out and place them in the freezer.

# PCR is done!

---



After PCR is done, you can either store your samples at  $-20^{\circ}\text{C}$  until you need to run a gel or to do DNA purification, or you can proceed immediately to DNA purification (~40 min total).



# Image and video credits

---

- Videos: PCR steps were demonstrated by Anya Wilkinson; recording and editing were done by Alina Avanesyan
- Photos: preparing and editing by Anya Wilkinson and Alina Avanesyan
- PCR protocol: prepared by Alina Avanesyan
- Videos were recorded and photos were taken in Dr. David Hawthorne's lab: 4172 Plant Science Building, Department of Entomology, University of Maryland, College Park, MD

## Acknowledgements

- We thank Dr. David Hawthorne (Department of Entomology, University of Maryland) for providing lab equipment and lab space for our DNA barcoding work; for providing lab space to take the photos and record the videos needed for developing this course; and for continuous support and encouragement!
- We also thank Anya Wilkinson for demonstrating the main PCR steps for this course