

## PCR protocol (plant DNA detection from leafhoppers)

1. Turn the thermocycler on (the black one on the left)
2. Take out from the freezer: PCR PreMix, primers (working solution; 4  $\mu$ M)
3. Prepare on the bench: ddH<sub>2</sub>O (in the right drawer), one 1.5 ml microcentrifuge tube.
4. Prepare PCR cocktail (#reactions = #samples+1; to cover pipetting error):

PCR reagent	1 reaction	5 reactions	10 reactions
PCR PreMix	10 $\mu$ l	50 $\mu$ l	100 $\mu$ l
Primer 1 (2 $\mu$ M)	2 $\mu$ l	10 $\mu$ l	20 $\mu$ l
Primer 2 (2 $\mu$ M)	2 $\mu$ l	10 $\mu$ l	20 $\mu$ l
ddH <sub>2</sub> O	5.2 $\mu$ l	26 $\mu$ l	52 $\mu$ l

5. Prepare 0.2 ml PCR strip tubes (in the drawer, next to the window)
6. Take DNA samples out from the fridge.
7. Place **19.2  $\mu$ l** of PCR cocktail in each tube.
8. Add **0.8  $\mu$ l** of a DNA sample in each tube. Mixed by pipetting.
9. Place strip tubes in the thermocycler (in any block: right or left)
10. Click on "RUN" -> "PROCEED". "Under "MAIN" find protocol ALINA1LH. Run it.

PCR settings:

- 94<sup>0</sup>C for 4 min
- $\times 35$  cycles:
  - 94<sup>0</sup>C for 30 sec
  - 57<sup>0</sup>C for 30 sec
  - 72<sup>0</sup>C for 30 sec
- 72<sup>0</sup>C for 2 min
- Hold at 4<sup>0</sup>C for 9 hours

11. To stop the protocol, press on "CANCEL". Turn the thermocycler off. Take the tubes out and place them in the freezer.

\* Annealing temperature and amount of DNA in a tube can be adjusted if needed.

\*\* After step 11 you can proceed immediately to gel electrophoresis if needed.