

Primer dilutions

Stock solution

Use #nmoles from the primer specification sheet.

For **100μM** : add (#nmoles × 10) μL of water

For **200μM** : add (#nmoles × 5) μL of water

Working solution

For **4 μM**:

$$C_1V_1 = C_2V_2$$

100μM (stock concentration) × **X** (stock volume needed) = 4 μM (working concentration) × **200 μL (working volume needed)**

$$\mathbf{X} = (4 \mu\text{M} \times 200 \mu\text{L}) / 100\mu\text{M}$$

$$\mathbf{X} = \mathbf{8 \mu\text{L}}$$

To prepare working solution:

$$\mathbf{8 \mu\text{L}} \text{ (stock solution)} + \mathbf{192 \mu\text{L}} \text{ (water)} = \mathbf{200 \mu\text{L}} \text{ (working volume needed)}$$